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## 15 Bibliographie

1. Kroll MH, Elinn RJ. Interference with Clinical Laboratory Analyses. *Clinical Chemistry*. 1994;40(11):1996–2005.
  2. Abramson N, Melton B. Leukocytosis: Basics of Clinical Assessment. *American Family Physician*. 2000 Nov 1;62(9):2053–60.
  3. Kurec A. Lipemia and hyperleukocytosis can lead to CBC errors. *Medical Laboratory Observer (MLO)*. 2016;48(3):44
  4. Riley LK, Rupert J. Evaluation of Patients with Leukocytosis. *American Family Physician*. 2015;92(11):1004–11.
  5. Apodaca MC, Wright AE, Riggins AM, Harris WP, Yeung RS, Yu L, et al. Characterization of a whole blood assay for quantifying myeloid-derived suppressor cells. *Journal for ImmunoTherapy of Cancer*. 2019;7(1).
  6. Constantino BT, Cogionis B. Nucleated RBCs—Significance in the Peripheral Blood Film. *Laboratory Medicine*. 2000;31(4):223–9.
  7. Buoro S, Vavassori M, Pipitone S, Benegiamo A, Lochis E, Fumagalli S, et al. Evaluation of nucleated red blood cell count by Sysmex XE-2100 in patients with thalassaemia or sickle cell anaemia and in neonates. *Blood Transfusion*. 2015;13(4):588–94.
  8. Booth F, Mead S v. Resistance to lysis of erythrocytes containing haemoglobin C-detected in a differential white cell counting system. *Journal of Clinical Pathology*. 1983; 36.
  9. Posteraro A, Gottfried EL. The diagnostic significance of a prolonged erythrocytic glycerol lysis time (GLT50). *American Journal of Clinical Pathology*. 1978;70(4):637–41.
  10. Genuardi E, Barbero D, Dogliotti I, Mantoan B, Drandi D, Gambella M, et al. Ficoll-hypaque separation vs whole blood lysis: Comparison of efficiency and impact on minimal residual disease analysis. *International Journal of Laboratory Hematology*. 2018;40(2):201–8.
  11. Dagur PK, McCoy JP. Collection, storage, and preparation of human blood cells. *Current Protocols in Cytometry*. 2015;2015:5.1.1-5.1.16.
  12. Buoro S, Manenti B, Seghezzi M. Which clinical significance has automatic detection of very low levels of nucleated red blood cells in the peripheral blood? *Annals of Translational Medicine*. 2016;4(11).
  13. Danise P, Amendola G, di Concilio R, Cillari E, Gioia M, di Palma A, et al. Nucleated red blood cells and soluble transferrin receptor in thalassemia syndromes: relationship with global and ineffective erythropoiesis. *Clinical Chemistry and Laboratory Medicine*. 2009;47(12):1539–42.
-

14. Lima M. Laboratory studies for paroxysmal nocturnal hemoglobinuria, with emphasis on flow cytometry. *Practical Laboratory Medicine*. 2020;20.
  15. Croom RD, McMillan CW, Sheldon GF, Orringer EP. Hereditary Spherocytosis Recent Experience and Current Concepts of Pathophysiology. *Annals of Surgery*. 1986;203(1):34–9.
  16. Wood B, Jevremovic D, Béné MC, Yan M, Jacobs P, Litwin V. Validation of Cell-based Fluorescence Assays: Practice Guidelines from the ICSH and ICCS – Part V – Assay Performance Criteria. *Cytometry Part B: Clinical Cytometry*. 2013;84(5):315–23.
  17. Farrell CJL, Carter AC. Serum indices: managing assay interference. Vol. 53, *Annals of Clinical Biochemistry*. SAGE Publications Ltd; 2016; 527–38.
  18. Yoo G, Kim J, Uh Y, Yoon KR, Park SD, Yoon KJ. Scoring system for detecting spurious hemolysis in anticoagulated blood specimens. *Annals of Laboratory Medicine*. 2015;35(3):341–7.
  19. Dimeski G. Interference Testing. *Clinical Biochemist Reviews*. 2008;29:43–8.
  20. Lippi G, Pavesi F, Benegiamo A, Pipitone S. What Do Hemolyzed Whole-Blood Specimens Look Like? Analysis with a CellaVision DM96 Automated Image Analysis System. *Journal of Laboratory Automation*. 2015;20(1):60–3.
  21. Weisbrot IM, Hollenberg LM. Platelet Counting Methods. *Laboratory Medicine*. 1980;11(5):307–12.
  22. Voulgaridou A, Kalfa TA. Autoimmune hemolytic anemia in the pediatric setting. *Journal of Clinical Medicine*. 2021;10(2):1–13.
  23. Lambert C, Yanikkaya Demirel G, Keller T, Preijers F, Psarra K, Schiemann M, et al. Flow Cytometric Analyses of Lymphocyte Markers in Immune Oncology: A Comprehensive Guidance for Validation Practice According to Laws and Standards. *Frontiers in Immunology*. 2020;11.
  24. García-González E, Aramendía M, Álvarez-Ballano D, Trincado P, Rello L. Serum sample containing endogenous antibodies interfering with multiple hormone immunoassays. *Laboratory strategies to detect interference*. *Practical Laboratory Medicine*. 2016;4:1–10.
  25. Rauch P, Zellmer A, CANDOR Bioscience GmbH Münster, Dankbar N, Institute of analytical chemistry university of Münster, Specht C, et al. Optimisation of assays: Interference in immunoassays recognize and avoid. *LABORWELT, Das BioTechnologie-Themenheft*. 2005;6(4):2–7.
  26. Selby C. Interference in immunoassay. Review Article *Annals of Clinical Biochemistry*. 1999;36:704–21.
  27. Tate J, Ward G. Interferences in Immunoassay. *Clinical Biochemist Reviews*. 2004;25:105–19.
  28. Luzzi VI, Scott MG, Gronowski AM. Negative Thyrotropin Assay Interference Associated with an IgGκ Paraprotein. *Clinical Chemistry*. 2003;49(4):709–10.
  29. Narayanan S. The Preanalytic Phase An Important Component of Laboratory Medicine. *American Journal of Clinical Pathology*. 2000;113:429–52.
  30. Cornbleet J. Spurious Results from Automated Hematology Cell Counters. *Laboratory Medicine*. 1983;14(8):509–14.
  31. Ghazal K, Brabant S, Prie D, Piketty ML. Hormone immunoassay interference: A 2021 update. *Annals of Laboratory Medicine*. 2021;42(1):3–23.
  32. Holm BE, Sandhu N, Tronstrøm J, Lydolph M, Trier NH, Houen G. Species cross-reactivity of rheumatoid factors and implications for immunoassays. *Scandinavian Journal of Clinical and Laboratory Investigation*. 2014;75(1):51–63.
  33. Dorn-Beineke A, Sack U. Quality control and validation in flow cytometry. *Journal of Laboratory Medicine - De Gruyter*. 2016;40(2):1–13.
  34. Renzi P, Ginns LC. Analysis of T cell subsets in normal adults Comparison of whole blood lysis technique to Ficoll-Hypaque separation by flow cytometry. *Journal of Immunological Methods*. 1987;98:53–6.
-

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Pacific Blue™ et Pacific Orange™ sont des marques commerciales de Life Technologies Corporation.